



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

KAPC A/B Antibody [#C10410]

Catalog Number: C10410

Concentration: 1mg/ml

Swiss-Prot No.: P17612

Other Names: EC 2.7.11.11, KAPCA, PKA C-alpha, PKA-alpha, PKACA, PRKACA, alpha-catalytic subunit, cAMP-dependent protein kinase, cAMP-dependent protein kinase, alpha-catalytic subunit, kinase PKA C-alpha

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human

KAPC A/B.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: KAPC A/B antibody detects endogenous levels of total KAPC A/B protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Maldonado F., Nucleic Acids Res. 16:8189-8190(1988).

Desseyn J.-L., Proc. Natl. Acad. Sci. U.S.A. 97:6433-6438(2000).

Ficarro S., J. Biol. Chem. 278:11579-11589(2003).

