



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

K6PL Antibody [#C10412]

Catalog Number: C10412

Concentration: 1mg/ml

Swiss-Prot No.: P17858

Other Names: 6-phosphofructokinase, liver type, EC 2.7.1.11, K6PL, PFK-B, PFKL, phosphofructo-1-kinase isozyme B, phosphofructokinase 1, phosphohexokinase

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human

K6PL.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: K6PL antibody detects endogenous levels of total K6PL protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Levanon D., DNA 8:733-743(1989).

Elson A., Genomics 7:47-56(1990).

Hattori M., Nature 405:311-319(2000).

