



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## EGR1/2 Antibody [#C10416]

**Catalog Number:** C10416

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P18146/P11161

**Other Names:** AT225, EGR-1, Early growth response protein 1, KROX-24, Krox-24 protein, NGFI-A, Nerve growth factor-induced protein A, Transcription factor ETR103, ZIF268, ZNF225, Zinc finger protein 225

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

EGR1/2.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** EGR1/2 antibody detects endogenous levels of total EGR1/2 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Suggs S.V., Nucleic Acids Res. 18:4283-4283(1990).

Shimizu N., J. Biochem. 111:272-277(1992).

Wright J.J., Science 248:588-591(1990).

