



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## ATF3 Antibody [#C10420]

**Catalog Number:** C10420

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P18847

**Other Names:** Cyclic AMP-dependent transcription factor  
ATF-3, Activating transcription factor 3, ATF3

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered  
saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl,  
0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against  
synthesized peptide derived from internal of human ATF3.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using  
epitope-specific immunogen.

**Specificity:** ATF3 antibody detects endogenous levels of  
total ATF3 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:5000

**References:**

Chen B.P.C., J. Biol. Chem. 269:15819-15826(1994).

Hai T., Genes Dev. 3:2083-2090(1989).

Hillman R.T., Genome Biol. 5:RESEARCH008.1-RESEARCH008.16(2004).

