



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## TSPAN8 Antibody [#C10423]

**Catalog Number:** C10423

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P19075

**Other Names:** Tetraspanin-8, Tspan-8, Transmembrane 4 superfamily member 3, Tumor-associated antigen CO-029, TSPAN8, TM4SF3

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human TSPAN8.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


**Specificity:** TSPAN8 antibody detects endogenous levels of total TSPAN8 protein.

**Reactivity:** Human


**Applications:** WB: 1:500~1:3000 IF: 1:100~1:500  
ELISA: 1:20000

**References:**

Szala S., Proc. Natl. Acad. Sci. U.S.A. 87:6833-6837(1990).  
The MGC Project Team, Genome Res. 14:2121-2127(2004).



Western blot analysis of extracts from 293 cells and COLO cells, using TSPAN8 antibody.



Immunofluorescence analysis of HeLa cells, using TSPAN8 antibody.