



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Elk1 Antibody [#C10425]

Catalog Number: C10425

Concentration: 1mg/ml

Swiss-Prot No.: P19419

Other Names: ELK1, ETS-domain protein ELK-1,
ETS-domain protein Elk-1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human Elk1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Elk1 antibody detects endogenous levels of total Elk1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Rao V.N., Science 244:66-70(1989).

Harindranath N., Gene 221:215-224(1998).

Yamauchi T., DNA Res. 6:21-27(1999).

