



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TFE3 Antibody [#C10430]

Catalog Number: C10430

Concentration: 1mg/ml

Swiss-Prot No.: P19532

Other Names: Transcription factor E3, TFE3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human TFE3.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TFE3 antibody detects endogenous levels of total TFE3 protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
ELISA: 1:20000

References:

Clark J., *Oncogene* 15:2233-2239(1997).


Sidhar S.K., *Hum. Mol. Genet.* 5:1333-1338(1996).

Beckmann H., *Genes Dev.* 4:167-179(1990).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using TFE3 antibody.



Western blot analysis of extracts from K562 cells and A549 cells, using TFE3 antibody.

For Research Use Only