



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## TBP Antibody [#C10435]

**Catalog Number:** C10435

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P20226

**Other Names:** TATA binding factor, TATA sequence-binding protein, TATA-box binding protein, TATA-box factor, TFIID, transcription initiation factor TFIID TBP subunit

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against

synthesized peptide derived from internal of human TBP.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** TBP antibody detects endogenous levels of total TBP protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:40000

**References:**

Peterson M.G., Science 248:1625-1630(1990).

Mungall A.J., Nature 425:805-811(2003).

Gangloff Y.G., Mol. Cell. Biol. 21:5109-5121(2001).

