



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TNF Receptor II Antibody [#C10436]

Catalog Number: C10436

Concentration: 1mg/ml

Swiss-Prot No.: P20333

Other Names: Tumor necrosis factor receptor superfamily member 1B [Precursor], Tumor necrosis factor receptor 2, TNF-R2, Tumor necrosis factor receptor type II, p75, p80 TNF-alpha receptor, CD120b antigen, Etanercept, Tumor necrosis factor receptor superfamily member 1b, membrane form, Tumor necrosis factor-binding protein 2, TNFRSF1B, TNFBR, TNFR2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from C-terminal of human TNF Receptor II.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TNF Receptor II antibody detects endogenous levels of total TNF Receptor II protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:20000

References:


Smith C.A., Science 248:1019-1023(1990).

Kohno T., Proc. Natl. Acad. Sci. U.S.A. 87:8331-8335(1990).

Beltinger C.P., Genomics 35:94-100(1996).



Western blot analysis of extracts from RAW264.7 cells, HT-29 cells and LOVO cells, using TNF Receptor II antibody.



Immunofluorescence analysis of HeLa cells, using TNF Receptor II antibody.

For Research Use Only