



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TNAP3 Antibody [#C10444]

Catalog Number: C10444

Concentration: 1mg/ml

Swiss-Prot No.: P21580

Other Names: A20, DNA binding protein A20, EC 3.-.-.-, TNAP3, TNFA1P2, tumor necrosis factor, alpha-induced protein 3, zinc finger protein A20

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human TNAP3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TNAP3 antibody detects endogenous levels of total TNAP3 protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:20000

References:


Opipari A.W. Jr., J. Biol. Chem. 265:14705-14708(1990).

Mungall A.J., Nature 425:805-811(2003).

Eliopoulos A.G., J. Virol. 73:1023-1035(1999).



Western blot analysis of extracts from RAW264.7 cells, using TNAP3 antibody.



Immunofluorescence analysis of HepG2 cells, using TNAP3 antibody.