



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

KAPCG Antibody [#C10451]

Catalog Number: C10451

Concentration: 1mg/ml

Swiss-Prot No.: P22612

Other Names: cAMP-dependent protein kinase, gamma-catalytic subunit, EC 2.7.11.11, PKA C-gamma, PRKACG

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human KAPCG.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: KAPCG antibody detects endogenous levels of total KAPCG protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Beebe S.J., Mol. Endocrinol. 4:465-475(1990).

Reinton N., Genomics 49:290-297(1998).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

