



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MYF6 Antibody [#C10457]

Catalog Number: C10457

Concentration: 1mg/ml

Swiss-Prot No.: P23409

Other Names: MRF4, Myf6, myogenic factor 6

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MYF6.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MYF6 antibody detects endogenous levels of total MYF6 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Braun T., EMBO J. 9:821-831(1990).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Kerst B., Hum. Mutat. 12:220-220(1998).

