



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

DGKA Antibody [#C10463]

Catalog Number: C10463

Concentration: 1mg/ml

Swiss-Prot No.: P23743

Other Names: 80 kDa diacylglycerol kinase, DAG kinase alpha, DAGK, DAGK1, DGK- alpha, DGKA, Diacylglycerol kinase, alpha, Diglyceride kinase, EC 2.7.1.107

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

DGKA.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: DGKA antibody detects endogenous levels of total DGKA protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Schaap D., FEBS Lett. 275:151-158(1990).

Hart T.C., Mamm. Genome 5:123-124(1994).

Hart T.C., Genomics 22:246-247(1994).

