



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ARBK1 Antibody [#C10475]

Catalog Number: C10475

Concentration: 1mg/ml

Swiss-Prot No.: P25098

Other Names: ADRBK1, ARBK1, BARK, BARK1, Beta-ARK-1, Beta-adrenergic receptor kinase 1, EC 2.7.11.15, G protein receptor kinase 2, G- protein coupled receptor kinase 2, kinase GRK2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

ARBK1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ARBK1 antibody detects endogenous levels of total ARBK1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Chuang T.T., J. Biol. Chem. 267:6886-6892(1992).
Penn R.B., J. Biol. Chem. 269:14924-14930(1994).
Fushman D., J. Biol. Chem. 273:2835-2843(1998).

