



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

PTX3 Antibody [#C10482]

Catalog Number: C10482

Concentration: 1mg/ml

Swiss-Prot No.: P26022

Other Names: Pentraxin-related protein PTX3
[Precursor], Pentaxin-related protein PTX3, Tumor
necrosis factor-inducible protein TSG-14, PTX3, TNFAIP5,
TSG14

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from N-terminal of human

PTX3.

Purification: The antibody was affinity-purified from
rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.


Specificity: PTX3 antibody detects endogenous levels of
total PTX3 protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:10000

References:

Breviario F., J. Biol. Chem. 267:22190-22197(1992).
Introna M., Eur. Heart J. 14:78-81(1993).
Lee G.W., J. Immunol. 150:1804-1812(1993).



Western blot analysis of extracts from 293 cells, using
PTX3 antibody.



Immunofluorescence analysis of HUVEC cells, using
PTX3 antibody.