



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ELAV2/4 Antibody [#C10486]

Catalog Number: C10486

Concentration: 1mg/ml

Swiss-Prot No.: Q12926/P26378

Other Names: ELAV-like protein 4, Paraneoplastic encephalomyelitis antigen HuD, Hu-antigen D, ELAVL4, HUD, PNEM

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ELAV2/4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ELAV2/4 antibody detects endogenous levels of total ELAV2/4 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:10000

References:


Szabo A., Cell 67:325-333(1991).

Behrends U., Int. J. Cancer 100:669-677(2002).

Gregory S.G., Nature 441:315-321(2006).



Western blot analysis of extracts from HUVEC cells and Jurkat cells, using ELAV2/4 antibody.



Immunofluorescence analysis of HepG2 cells, using ELAV2/4 antibody.