



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

APEX1 Antibody [#C10493]

Catalog Number: C10493

Concentration: 1mg/ml

Swiss-Prot No.: P27695

Other Names: AP endonuclease 1, APE, APE1, APEN, APEX, APEX nuclease, APEX1, APX, DNA-(apurinic or apyrimidinic site) lyase, EC 4.2.99.18, HAP1, REF-1 protein, REF1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

APEX1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: APEX1 antibody detects endogenous levels of total APEX1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Robson C.N., Nucleic Acids Res. 20:4417-4421(1992).

Xanthoudakis S., EMBO J. 11:3323-3335(1992).

Zhao B., Nucleic Acids Res. 20:4097-4098(1992).

