



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GPR126 Antibody [#G290]

Catalog Number: G290

Concentration: 1mg/ml

Swiss-Prot No.: Q86SQ4

Other Names: GP126; GPR126

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GPR126.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GPR126 antibody detects endogenous levels of total GPR126 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Ota T., Nat. Genet. 36:40-45(2004).

Fredriksson R., Biochem. Biophys. Res. Commun. 301:725-734(2003).

Kristiansen T.Z., Mol. Cell. Proteomics 3:715-728(2004).

