



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** +86 519 8805 0026

## GPR144 Antibody [#G298]

**Catalog Number:** G298

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q7Z7M1

**Other Names:** Probable G-protein coupled receptor 144;  
G-protein coupled receptor PGR24; GPR144; PGR24

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human GPR144.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** GPR144 antibody detects endogenous levels of total GPR144 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:1000

**References:**

Bjarnadottir T.K., Genomics 84:23-33(2004).

Humphray S.J., Nature 429:369-374(2004).

Vassilatis D.K., Proc. Natl. Acad. Sci. U.S.A. 100:4903-4908(2003).

