



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GPR150 Antibody [#G301]

Catalog Number: G301

Concentration: 1mg/ml

Swiss-Prot No.: Q8NGU9

Other Names: Probable G-protein coupled receptor 150

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GPR150.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GPR150 antibody detects endogenous levels of total GPR150 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 92%, Positives = 92%)

Applications: IF: 1:100~1:500 ELISA: 1:1000

References:

Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.

