



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GPR152 Antibody [#G303]

Catalog Number: G303

Concentration: 1mg/ml

Swiss-Prot No.: Q8TDT2

Other Names: Probable G-protein coupled receptor 152;
G-protein coupled receptor PGR5; GPR152; PGR

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GPR152.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GPR152 antibody detects endogenous levels of total GPR152 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:1000

References:

Bonner T.I., Submitted (MAR-2004) to the EMBL/GenBank/DDBJ databases.

Takeda S., FEBS Lett. 520:97-101(2002).

Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.

