



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GPR153 Antibody [#G304]

Catalog Number: G304

Concentration: 1mg/ml

Swiss-Prot No.: Q6NV75

Other Names: DKFZp762B2210; G-protein coupled receptor 153; G-protein coupled receptor PGR1; GP153; PGR1v

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GPR153.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GPR153 antibody detects endogenous levels of total GPR153 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:5000

References:


Gregory S.G., Nature 441:315-321(2006).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Vassilatis D.K., Proc. Natl. Acad. Sci. U.S.A. 100:4903-4908(2003).



Western blot analysis of extracts from K562 cells, using GPR153 antibody.



Immunofluorescence analysis of A549 cells, using GPR153 antibody.