



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GPR157 Antibody [#G306]

Catalog Number: G306

Concentration: 1mg/ml

Swiss-Prot No.: Q5UAW9

Other Names: Probable G-protein coupled receptor 157

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GPR157.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GPR157 antibody detects endogenous levels of total GPR157 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:5000

References:


Bonner T.I., Submitted (SEP-2004) to the EMBL/GenBank/DDBJ databases.

Ota T., Nat. Genet. 36:40-45(2004).

Gregory S.G., Nature 441:315-321(2006).



Western blot analysis of extracts from NIH-3T3 cells, using GPR157 antibody.



Immunofluorescence analysis of HUVEC cells, using GPR157 antibody.