



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GPR175 Antibody [#G313]

Catalog Number: G313

Concentration: 1mg/ml

Swiss-Prot No.: Q86W33

Other Names: Integral membrane protein GPR175

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GPR175.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Specificity: GPR175 antibody detects endogenous levels of total GPR175 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:1000

References:

Saito T., Saito T., Submitted (JAN-2000) to the EMBL/GenBank/DDBJ databases.

Wan D., Proc. Natl. Acad. Sci. U.S.A. 101:15724-15729(2004).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

