



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ERCC5 Antibody [#C10503]

Catalog Number: C10503

Concentration: 1mg/ml

Swiss-Prot No.: P28715

Other Names: DNA excision repair protein ERCC-5, DNA-repair protein complementing XP-G cells, ERCC5, XPGC, Xeroderma pigmentosum group G complementing protein

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

ERCC5.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ERCC5 antibody detects endogenous levels of total ERCC5 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Scherly D., Nature 363:182-185(1993).

Shiomi T., Mutat. Res. 314:167-175(1994).

Macinnes M.A., Mol. Cell. Biol. 13:6393-6402(1993).

