



Web: [www.anbobio.com](http://www.anbobio.com)

[www.dijibio.com](http://www.dijibio.com)

Order: [order@dijibio.com](mailto:order@dijibio.com)

Tel: (86) 519 8805 0026

## EPHA3 Antibody [#C10508]

**Catalog Number:** C10508

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P29320

**Other Names:** EPA3; ETK; ETK1; Ephrin type-A receptor 3 precursor; HEK; HEK4; MEK4; REK4; TYRO-4; TYRO4; Tyrosine-protein kinase receptor ETK1; Tyrosine-protein kinase receptor REK4; kinase EphA3

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human EPHA3.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** EPHA3 antibody detects endogenous levels of total EPHA3 protein.

**Reactivity:** Human (Identities = 100%, Positives = 100%)

Mouse (Identities = 100%, Positives = 100%)

Rat (Identities = 100%, Positives = 100%)

**Applications:** WB: 1:500~3000 IHC: 1:50~1:100

ELISA: 1:40000

**References:**

Coulie P.G., Cancer Res. 60:4855-4863(2000).

Comb M.J., Cell 131:1190-1203(2007)

Brenner S.E., Genome Biol. 5:R8.1-R8.16(2004).

