



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

EPHB2 Antibody [#C10510]

Catalog Number: C10510

Concentration: 1mg/ml

Swiss-Prot No.: P29323

Other Names: CEK5, DRT, EC 2.7.10.1, EPH-3, EPHB2, EPTH3, ERK, Ephrin type-B receptor 2 precursor, HEK5, NUK, Receptor protein-tyrosine kinase HEK5, SEK-3, SEK3, kinase EphB2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human

EPHB2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: EPHB2 antibody detects endogenous levels of total EPHB2 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Kiyokawa E., Cancer Res. 54:3645-3650(1994).

Ikegaki N., Hum. Mol. Genet. 4:2033-2045(1995).

Tang X.X., Oncogene 17:521-526(1998).

