



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

LTK Antibody [#C10512]

Catalog Number: C10512

Concentration: 1mg/ml

Swiss-Prot No.: P29376

Other Names: EC 2.7.1.112, EC 2.7.10.1, Leukocyte tyrosine kinase receptor precursor, Protein tyrosine kinase-1, TYK1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human LTK.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: LTK antibody detects endogenous levels of total LTK protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Toyoshima H., Proc. Natl. Acad. Sci. U.S.A. 90:5404-5408(1993).

Krolewski J.J., EMBO J. 10:2911-2919(1991).

Maru Y., Oncogene Res. 5:199-204(1990).

