



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

KCY Antibody [#C10514]

Catalog Number: C10514

Concentration: 1mg/ml

Swiss-Prot No.: P30085

Other Names: CMK, EC 2.7.4.14, KCY, RP11-511I2.1, UCK, UMK, UMP-CMP kinase, UMP-CMPK, UMP/CMP kinase, UMP/CMPK, UMPK, cytidine monophosphate kinase, cytidylate kinase, deoxycytidylate kinase, uridine monophosphate kinase, uridine monophosphate/cytidine monophosphate kinase

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human KCY.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: KCY antibody detects endogenous levels of total KCY protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Van Rompay A.R., Mol. Pharmacol. 56:562-569(1999).
Pearman A.T., Life Sci. 69:2361-2370(2001).
Liou J.-Y., Cancer Res. 62:1624-1631(2002).

