



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

KAP1 Antibody [#C10527]

Catalog Number: C10527

Concentration: 1mg/ml

Swiss-Prot No.: P31321

Other Names: KAP1, PRKAR1B, cAMP-dependent protein kinase type I-beta regulatory subunit

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human KAP1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: KAP1 antibody detects endogenous levels of total KAP1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Solberg R., Biochem. Biophys. Res. Commun. 176:166-172(1991).

Solberg R., Exp. Cell Res. 214:595-605(1994).

Zhan X., Anal. Biochem. 354:279-289(2006).

