



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

GLPK Antibody [#C10533]

Catalog Number: C10533

Concentration: 1mg/ml

Swiss-Prot No.: P32189

Other Names: Glycerol kinase, EC 2.7.1.30, ATP:glycerol 3-phosphotransferase, Glycerokinase, GKI

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human GLPK.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GLPK antibody detects endogenous levels of total GLPK protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:20000

References:


Guo W., Nat. Genet. 4:367-372(1993).

Sargent C.A., Hum. Mol. Genet. 3:1317-1324(1994).

Walker A.P., Hum. Mol. Genet. 2:107-114(1993).



Western blot analysis of extracts from 293 cells and Jurkat cells, using GLPK antibody.



Immunofluorescence analysis of HeLa cells, using GLPK antibody.