



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

FOXN2 Antibody [#C10537]

Catalog Number: C10537

Concentration: 1mg/ml

Swiss-Prot No.: P32314

Other Names: Forkhead box protein N2, Human T-cell leukemia virus enhancer factor, FOXN2, HTLF

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human FOXN2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FOXN2 antibody detects endogenous levels of total FOXN2 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Li C., Submitted (APR-1996) to the EMBL/GenBank/DDBJ databases.

Li C., Genomics 13:658-664(1992).

