



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ELF1 Antibody [#C10538]

Catalog Number: C10538

Concentration: 1mg/ml

Swiss-Prot No.: P32519

Other Names: E74-like factor 1, ELF1, ETS-related transcription factor Elf-1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human ELF1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ELF1 antibody detects endogenous levels of total ELF1 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:1000

References:


Leiden J.M., J. Virol. 66:5890-5897(1992).

Thompson C.B., Mol. Cell. Biol. 12:1043-1053(1992).

Wang C.Y., Science 260:1330-1335(1993).



Western blot analysis of extracts from 293 cells and HuvEc cells, using ELF1 antibody.



Immunofluorescence analysis of A549 cells, using ELF1 antibody.