



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GPRC5C Antibody [#G319]

Catalog Number: G319

Concentration: 1mg/ml

Swiss-Prot No.: Q9NQ84

Other Names: GC5C; GPC5C; RAIG-3; RAIG3, G protein-coupled receptor family C group 5 member C precursor; retinoic acid induced gene 3 protein

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GPRC5C.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GPRC5C antibody detects endogenous levels of total GPRC5C protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Robbins M.J., Genomics 67:8-18(2000).

Braeuner-Osborne H., Biochim. Biophys. Acta 1518:237-248(2001).

Matsuoka S., Science 316:1160-1166(2007).

