



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GPC6A Antibody [#G321]

Catalog Number: G321

Concentration: 1mg/ml

Swiss-Prot No.: Q5T6X5

Other Names: G protein-coupled receptor, family C, group 6, member A; G-protein coupled receptor 33; GPC6A; GPCR33

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GPC6A.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GPC6A antibody detects endogenous levels of total GPC6A protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:1000

References:


Wellendorph P., Gene 335:37-46(2004).

Mungall A.J., Nature 425:805-811(2003).

Takeda S., FEBS Lett. 520:97-101(2002).



Western blot analysis of extracts from Jurkat cells, using GPC6A antibody.



Immunofluorescence analysis of MCF-7 cells, using GPC6A antibody.