



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GPR17 Antibody [#G325]

Catalog Number: G325

Concentration: 1mg/ml

Swiss-Prot No.: Q13304

Other Names: Uracil nucleotide/cysteinyl leukotriene receptor; UDP/CysLT receptor; G-protein coupled receptor 17; P2Y-like receptor; R12; GPR17

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mMNaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

GPR17.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GPR17 antibody detects endogenous levels of total GPR17 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Raport C.J., J. Leukoc. Biol. 59:18-23(1996).

Blaesius R.H., J. Neurochem. 70:1357-1365(1998).

Ciana P., EMBO J. 25:4615-4627(2006).

