



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GPR19 Antibody [#G327]

Catalog Number: G327

Concentration: 1mg/ml

Swiss-Prot No.: Q15760

Other Names: Probable G-protein coupled receptor 19;
GPR-NGA; GPR19

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GPR19.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GPR19 antibody detects endogenous levels of total GPR19 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Bonner T.I., Submitted (MAY-1996) to the EMBL/GenBank/DDBJ databases.

O'Dowd B.F., FEBS Lett. 394:325-329(1996).

