



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GPR25 Antibody [#G331]

Catalog Number: G331

Concentration: 1mg/ml

Swiss-Prot No.: O00155

Other Names: Probable G-protein coupled receptor 25

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GPR25.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GPR25 antibody detects endogenous levels of total GPR25 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Jung B.P., Biochem. Biophys. Res. Commun. 230:69-72(1997).

Jung B.P., Submitted (NOV-2003) to the EMBL/GenBank/DDBJ databases.

