



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GPR31 Antibody [#G334]

Catalog Number: G334

Concentration: 1mg/ml

Swiss-Prot No.: O00270

Other Names: Probable G-protein coupled receptor 31;
GPR31

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

GPR31.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GPR31 antibody detects endogenous levels of total GPR31 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Zingoni A., Genomics 42:519-523(1997).

Mungall A.J., Nature 425:805-811(2003).

