



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GPR34 Antibody [#G337]

Catalog Number: G337

Concentration: 1mg/ml

Swiss-Prot No.: Q9UPC5

Other Names: G protein-coupled receptor 34

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GPR34.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GPR34 antibody detects endogenous levels of total GPR34 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:10000

References:

Schoneberg T., *Biochim. Biophys. Acta* 1446:57-70(1999).

Marchese A., *Genomics* 56:12-21(1999).

Jacobi F.K., *Hum. Genet.* 107:89-91(2000).



Western blot analysis of extracts from K562 cells, HT-29 cells and COS-7 cells, using GPR34 antibody.



Immunofluorescence analysis of LOVO cells, using GPR34 antibody.