



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** +86 519 8805 0026

## GPR37 Antibody [#G339]

**Catalog Number:** G339

**Concentration:** 1mg/ml

**Swiss-Prot No.:** O15354

**Other Names:** Probable G-protein coupled receptor 37;

Endothelin B receptor-like protein 1; ETBR-LP-1;

Parkin-associated endothelin receptor-like receptor;

PAELR; GPR37

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human GPR37.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** GPR37 antibody detects endogenous levels

of total GPR37 protein.

**Reactivity:** Human (Identities = 100%, Positives = 100%);

Mouse (Identities = 87%, Positives = 100%);

Rat (Identities = 93%, Positives = 100%)

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Marazziti D., Genomics 45:68-77(1997).

Donohue P.J., Brain Res. Mol. Brain Res. 54:152-160(1998).

Zeng Z., Biochem. Biophys. Res. Commun. 233:559-567(1997).

