



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CIITA Antibody [#C10541]

Catalog Number: C10541

Concentration: 1mg/ml

Swiss-Prot No.: P33076

Other Names: MHC class II transactivator, CIITA, CIITA, MHC2TA

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CIITA.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CIITA antibody detects endogenous levels of total CIITA protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Steimle V., Cell 75:135-146(1993).

Al-Kandari W., J. Mol. Biol. 369:1175-1187(2007).

Al-Kandari W., Mol. Immunol. 44:311-321(2007).

Croager E.J., Biochim. Biophys. Acta 1353:231-235(1997).

