



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

Catenin- β 1 Antibody [#C10546]

Catalog Number: C10546

Concentration: 1mg/ml

Swiss-Prot No.: P35222

Other Names: CATNB, CTNB1, CTNNB, catenin, beta

Storage/Stability: Store at -20°C /1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human Catenin- β 1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Catenin- β 1 antibody detects endogenous levels of total Catenin- β 1 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
ELISA: 1:10000

References:

Huelsken J., J. Cell Biol. 127:2061-2069(1994).


Kim J.-S., Cancer Res. 62:2744-2748(2002).

Kas K., Nat. Genet. 15:170-174(1997).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue using Catenin- β 1 antibody.



Western blot analysis of extracts from HeLa cells, using Catenin- β 1 antibody.