



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## Recoverin Antibody [#C10549]

**Catalog Number:** C10549

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P35243

**Other Names:** Recoverin, Cancer-associated retinopathy protein, Protein CAR, RCVRN, RCV1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human Recoverin.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** Recoverin antibody detects endogenous levels of total Recoverin protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Murakami A., Biochem. Biophys. Res. Commun. 187:234-244(1992).

Wiechmann A.F., Exp. Eye Res. 56:463-470(1993).

Thirkill C.E., Invest. Ophthalmol. Vis. Sci. 33:2768-2772(1992).

