



Web: [www.anbobio.com](http://www.anbobio.com)

[www.dijibio.com](http://www.dijibio.com)

Order: [order@dijibio.com](mailto:order@dijibio.com)

Tel: (86) 519 8805 0026

## GRK3 Antibody [#C10560]

**Catalog Number:** C10560

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P35626

**Other Names:** ADRBK2, ARBK2, BARK2, EC 2.7.11.15, G-protein coupled receptor kinase 3, beta-ARK-2, beta-adrenergic receptor kinase 2

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

GRK3.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** GRK3 antibody detects endogenous levels of total GRK3 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Parruti G., Biochem. Biophys. Res. Commun. 190:475-481(1993).

Dunham I., Nature 402:489-495(1999).

