



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

SOX6 Antibody [#C10565]

Catalog Number: C10565

Concentration: 1mg/ml

Swiss-Prot No.: P35712

Other Names: Transcription factor SOX-6, SRY-BOX 6

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human SOX6.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: SOX6 antibody detects endogenous levels of total SOX6 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Cohen-Barak O., Gene 265:157-164(2001).

Wiemann S., Genome Res. 11:422-435(2001).

Denny P., Nucleic Acids Res. 20:2887-2887(1992).

