



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## TCF7 Antibody [#C10571]

**Catalog Number:** C10571

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P36402

**Other Names:** Transcription factor 7, T-cell-specific transcription factor 1, TCF-1, T-cell factor 1, TCF7, TCF1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from N-terminal of human TCF7.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** TCF7 antibody detects endogenous levels of total TCF7 protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

van de Wetering M., EMBO J. 10:123-132(1991).

Mayer K., Biochim. Biophys. Acta 1263:169-172(1995).

Roose J., Science 285:1923-1926(1999).

