



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

CXCR3 Antibody [#G063]

Catalog Number: G063

Concentration: 1mg/ml

Swiss-Prot No.: P49682

Other Names: C-X-C chemokine receptor type 3; CD183; CD183 antigen; CKR-L2; CMKAR3; CXC-R3; CXCR-3; G protein-coupled receptor 9; GPR9; interferon-inducible protein 10 receptor; IP-10 receptor; IP10; IP10 receptor; IP10-R; Mig receptor; Mig-R; MigR

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from internal of human CXCR3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CXCR3 antibody detects endogenous levels of total CXCR3 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Loetscher M., J. Exp. Med. 184:963-969(1996).

Lasagni L., J. Exp. Med. 197:1537-1549(2003).

Marchese A., Genomics 29:335-344(1995).



For Research Use Only