



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

CCKBR Antibody [#G069]

Catalog Number: G069

Concentration: 1mg/ml

Swiss-Prot No.: P32239

Other Names: CCK-B receptor; CCK-BR; CCKRB;
Cholecystokinin-2 receptor; GASR;

Gastrin/cholecystokinin type B receptor

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CCKBR.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CCKBR antibody detects endogenous levels of total CCKBR protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Pisegna J.R., Biochem. Biophys. Res. Commun. 189:296-303(1992).

Lee Y.-M., J. Biol. Chem. 268:8164-8169(1993).

Ito M., J. Biol. Chem. 268:18300-18305(1993).

