



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

CLPT1 Antibody [#G070]

Catalog Number: G070

Concentration: 1mg/ml

Swiss-Prot No.: O96005

Other Names: Cleft lip and palate transmembrane protein 1; CLPTM1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CLPT1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CLPT1 antibody detects endogenous levels of total CLPT1 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Yoshiura K., Genomics 54:231-240(1998).

Totoki Y., Submitted (APR-2005) to the EMBL/GenBank/DDBJ databases.

Lewandrowski U., Mol. Cell. Proteomics 5:226-233(2006).

